



FACULTY OF PHARMACY  
IN HRADEC KRÁLOVÉ  
Charles University

25. 02. / 04. 03. / 11. 3. 2022


# Selected Methods of Pharmaceutical Technology

## Dissolution

Mgr. Jana Brokešová (brokesoj@faf.cuni.cz)

= determination of a **mass of released drug** from a dosage form in predetermined liquid and at predetermined time

Principle – schema of flow-through cell (Apparatus 4)





- ONLINE X OFFLINE analysis

Automatic sample analysis (during dissolution testing)


Manual sample analysis (after dissolution testing is complete)

- OPEN X CLOSE system



## Pharmacokinetics

The principles of ADME




- Why we are dealing with dissolution?
- dissolution x solubility

## Why?

- > dissolution profiles allow us to assume bioavailability

**Drug dissolving** – necessary for the absorption (and effect) of the drug in the organism

*FT I lecture: Biopharmaceutical aspects of dosage forms (Dr. Vraníková)*



## What is it good for?


- drug development and control of drug quality
  - study of increasing a solubility of API (e.g.: salt formation, surfactants, complexation, nanonization, amorphization, solid dispersions...)
  - bioequivalence study of the original and generic preparations
  - quality control (batch uniformity monitoring)

## Revision

*FT I lecture: Dissolution testing (Dr. Šnejdrová)*

- Four types of dissolution apparatus (Ph. Eur.)
- Apparent x Intrinsic dissolution
- Dissolution media (HCl, phosphate, etc.... purified water is used only if the release of a drug is not pH dependent)
- Temperature of medium during dissolution testing of solid dosage forms ( $37 \pm 0,5 \text{ }^\circ\text{C}$ )
- Sink conditions
  - during today's laboratory practice: flow-through cell apparatus

# Flow-through cell dissolution (Apparatus 4)



- a suitable cell is selected according to the type of dosage form
- medium (composition, volume, temperature)
- flow rate of medium
- collecting samples in specific time intervals x online analysis (apparatus with special flow-through cuvette)




*different types of cells and industry dissolution apparatus – Sotax seminar 8. 4. 2022*

## Advantages of flow-through cell dissolution apparatus


- Most versatile (tablets, powders, semisolids, contact lenses, stents, ...) (no previous compression of powders required)
- More accurate evaluation of the dissolution test (more sensitive method – captures smaller changes, better simulations of intraluminal hydrodynamics)
- Suitable if pH change is necessary (change of medium during the experiment) – evaluation of the dosage forms with controlled release
- Apparatus can be used to measure apparent dissolution – this method is mainly used to determine the apparent dissolution rate of pure solid substances (cell for powders) (Ph. Eur. 2.9.43)
- In combination with other tests / analysis – broaden the knowledge about the dosage form (research)

## Procedure – step by step

- choose a suitable cell – 2 sizes of cell for solid dosage forms, necessary to know the volume of the cell for the following calculations (ml)
- one ruby ball (5 mm) + glass balls (1 mm) – to ensure laminar flow in the cell (without glass balls – turbulent flow)
- cell composition – according to the type of dosage form






- placed the sample
- medium preparation - deaeration (UZ / heating+filtration)
- medium flow rate – 4, 8, 16, ... (ml/min)
- open system (suitable if the solubility of the drug is very low) / close system; offline / online analysis





- What the actual flow-through cell dissolution apparatus looks like:




Connection of the tubes (flow-through cell dissolution, close system):


medium -> pump > flow-through cuvette -> cell containing sample-> medium

cell containing tablet before the dissolution test

cell containing tablet during the dissolution testing



- What the actual basket and paddle apparatus look like:



- you can see off-line system, there is no connection with spectrophotometer, the measurement of the samples is manual, after the dissolution test is over

- if the tablets are too „light“ and they float on the surface, it is appropriate to use sinkers



*Basket and paddle apparatus was taught during the practical training of Pharmaceutical Technology*


# Evaluation

- preparation of standard solution sequence of known concentration -> calibration curve (c, A)
- measurement of collected samples absorbance (A)
- calculations:
  - concentration of samples (c) according to the calibration curve
  - mass of released drug at individual intervals ( $m_{int}$ ) according to this formula

$$m_{int} (mg) = c * V$$

(„V“ is volume of the medium)

- relative mass of released drug % ( $m_{rel}$ ) equals  $m_{int}$  divided by mass of the drug in the tested sample



## Graphs

open system - > non-cumulative curve

close system - > cumulative curve

## Ph. Eur.

complete release = at least 80 % API, till 45 min (conventional-release)

*FT I lecture: Dissolution testing (Dr. Šnejdrová)*



## Pharmacokinetics

- mathematical models are used for more detailed study


# Example of results evaluation (tablets containing paracetamol, measured by two methods – Paddle apparatus, Flow-through cell apparatus)

Time (min.)	A	$m_{rel}$ PAR (%)
0	0.000	0.00
2	0.000	2.72
4	0.010	3.78
6	0.024	5.33
8	0.059	9.18
10	0.100	13.72
12	0.129	16.91
14	0.157	20.06
16	0.196	24.34
18	0.226	27.64
20	0.278	33.46
22	0.311	37.11
24	0.329	39.10
26	0.353	41.73
28	0.383	45.09
30	0.414	48.46
32	0.444	51.83
34	0.475	55.20
36	0.505	58.56
38	0.536	61.93
40	0.566	65.30
42	0.596	68.66
44	0.627	72.03
46	0.657	75.40
48	0.688	78.77
50	0.718	82.13
52	0.749	85.50
54	0.779	88.87
56	0.809	92.24
58	0.840	95.60
60	0.870	98.97

Time (min.)	A	$m_{rel}$ PAR (%)
0	0.000	0.00
2	0.134	17.50
4	0.181	22.66
6	0.212	26.13
8	0.255	30.85
10	0.299	35.76
15	0.412	48.30
20	0.584	67.27
25	0.728	83.23
30	0.790	90.09
40	0.860	97.84
50	0.870	98.95
60	0.875	99.50



# Comparison of both methods – Paddle apparatus, Flow-through cell apparatus



Discussion:

What influences dissolution?

- dosage form, excipients
- chosen type of dissolution apparatus (different dissolution profiles)
- conditions:
  - stirrer speed
  - medium flow rate
  - medium volume
  - type of medium
  - pH and ionic strength of medium
- + properties of API (crystal (amorphous) structure, particle size, surface properties...)

*Physical Principles of Dosage Forms lecture: Solubility and Dissolution (Dr. Smékalová)*



FACULTY OF PHARMACY  
IN HRADEC KRÁLOVÉ  
Charles University

Thank you for your attention.

## Sources:

- Ph. Eur.
- faf.cuni.cz - Study materials
- Fotaki, N.: Flow-Through Cell Apparatus (USP Apparatus 4): Operation and Features, Dissolution Technologies, 2011  
[http://dissolutiontech.com/DTresour/201111Articles/DT201111\\_A06.pdf](http://dissolutiontech.com/DTresour/201111Articles/DT201111_A06.pdf)
- Cupera, J., Lansky, P., Sklubalova, Z.: Sampling times influence the estimate of parameters in the Weibull dissolution model, Eur. J. Pharm. Sci., 2015  
<https://www.ncbi.nlm.nih.gov/pubmed/26215461>
- <http://www.collegeofpharmacy.com/callearningcentre/demonstrationCALPackages/DissolutionDemo/1Introduction/page-7.html>
- <https://accesspharmacy.mhmedical.com/content.aspx?bookid=513&sectionid=41488032>
- <http://edelweisspublications.com/keyword/33/355/Pharmacokinetics>