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Login

BRENDA

GLOBAL CORE BIODATA RESOURCE **BRENDA has been selected as Global Core Biodata Resource** [More...](#)

Search term **Search BRENDA**

Enzyme, Ligand contains

[add search field](#) [delete search field](#)

- Text-based queries**
 - Full-text Search
 - Advanced Search
 - DRENDA
 - MeSH Ontology
- Structure-based queries**
 - Ligand Structure Search
 - Metabolic Pathways
 - Enzyme Structures
- Explorer**
 - Enzyme Classification
 - TaxTree
 - Protein folding: CATH, SCOPe Ontologies
- Visualization**
 - Enzyme Word Maps
 - Organism Word Maps
 - Genomes
 - Functional Parameter Statistics
 - Metabolic Pathways
- Prediction**
 - Membrane Helices
 - Localization Prediction
 - EnzymeDetector
- Supporting & External**
 - BRENDA
 - Tissue Ontology
 - Biochemical Reactions
 - MetaboMAPS
- News**
 - Introducing: Quick Feedback!** - September 23, 2022
Feedback from our users is very important to us! We want to make it as easy as possi...
 - New user system** - September 19, 2022
We have enabled a new login and user system. This should resolve problems with un...
 - BRENDA json download** - September 12, 2022
BRENDA is now available to download in json format!
The download can be found [here](#) and a documentation of the format is available ([here](#)...

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
2.6.1.1: aspartate transaminase

This is an abbreviated version!
For detailed information about aspartate transaminase, go to the full flat file.

Word Map on EC 2.6.1.1 

Reaction

Synonyms

More, Pat, glutamic-oxaloacetic transaminase, glutamic-aspartic transaminase, transaminase A, AAT, AspT, AspAT, aminotransferase, aspartate, 2-oxoglutarate-glutamate aminotransferase, aspartate alpha-ketoglutarate transaminase, aspartate aminotransferase, aspartate-2-oxoglutarate transaminase, aspartic acid aminotransferase, aspartic aminotransferase, aspartyl aminotransferase, AST, glutamate-oxalacetate aminotransferase, glutamate-oxalate transaminase, glutamic-aspartic aminotransferase, glutamic-oxalacetic transaminase, glutamic oxalic transaminase, GOT (enzyme), GOT, L-aspartate transaminase, L-aspartate-alpha-ketoglutarate transaminase, L-aspartate-2-ketoglutarate aminotransferase, L-aspartate-2-oxoglutarate aminotransferase, L-aspartate-2-oxoglutarate-transaminase, L-aspartate-2-oxoglutarate-transaminase, L-aspartic aminotransferase, oxaloacetate-aspartate aminotransferase, oxaloacetate transferase, aspartate:2-oxoglutarate aminotransferase, glutamate oxaloacetate transaminase, aspartate:2-oxoglutarate aminotransferase, L-aspartate aminotransferase, protein TT0402, AsAT, AAT3, plastid aspartate aminotransferase, aspartate, 2-oxoglutarate aminotransferase, prephenate aminotransferase, AATase, Sar2028, aspartate/tyrosine/phenylalanine pyridoxal-5'-phosphate-dependent aminotransferase, KAT IV, mitAAT, L-aspartate:2-oxoglutarate aminotransferase, aspartate transaminase, glutamic oxaloacetic transaminase, AspATs, SsAspAT, aspB, C-S lyase, aatA, CAA1, aspartate/(R)-cysteate:2-oxoglutarate aminotransferase, aoa/coa, PAspAT, L-AspAT, L-aspartate:2-oxoglutarate aminotransferase, EcAspAT, aatB3, AST-Bb, GOT1L1, alr1039, all2340, alr2765, alr4853, Cgl0240, aspartate aminotransferase A, Tb11.02.2740, GL50803_91056, GOT1, GOT2, AAT-2, AAT-3, Rv3722c, cAST, aspartate AT, Asp AT, AtPAT, bifunctional aspartate aminotransferase and glutamate/aspartate-prephenate aminotransferase, class I beta AAT, aspartate aminotransferase 1, glutamate oxaloacetate transaminase 1, PT-AAT, aspC, glutamate-oxaloacetate transaminase 1

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AST IFCC Multi-Purpose (MPR) Liquid Reagent

KIT SPECIFICATIONS:

Cat. No	Quantity	Reagent	Storage
GL733AS	10 x 15 ml	AST IFCC - 1	2-8°C
	2 x 15 ml	AST IFCC - 2	
GL743AS	5 x 50 ml	AST IFCC - 1	2-8°C
	1 x 50 ml	AST IFCC - 2	
GL753AS	5 x 100 ml	AST IFCC - 1	2-8°C
	1 x 100 ml	AST IFCC - 2	

INTENDED USE:

In Vitro Diagnostic reagent pack for the quantitative determination of Aspartate Aminotransferase (AST) in serum and plasma, based upon the IFCC recommendations, on automated and semi-automated analyzers.

SUMMARY AND EXPLANATION:

Aspartate aminotransferase (AST) belongs to a group of enzymes that catalyse the interconversion of amino acids to 2-oxo-acids. The enzyme is largely present in the cytoplasm of cells and highest concentrations are from heart, liver and skeletal muscle tissues. The main clinical uses are in the diagnosis and treatment of myocardial infarction and hepatic diseases. After myocardial infarction, AST level is elevated 5 to 8 hours after the onset of chest pain. It remains elevated for 4 days after which it returns to normal. Five to ten-fold elevation of AST are observed in primary metastatic carcinoma of the liver.

Pyridoxyl-5-Phosphate
Pyridoxyl-5-Phosphate functions as a co-enzyme in the aminotransferase reaction. While most patient samples contain endogenous pyridoxyl phosphate the IFCC recommendation specify the addition of the co-enzyme to supplement pyridoxyl-5-phosphate deficient samples. The pyridoxyl-5-phosphate is available as a powder vial from Glenbio Ltd. Cat. No. GL806PP (12 x 50 ml) and GL816PP (6 x 66 ml).

PRINCIPLE OF THE TEST: 1,2

The sample is pre-incubated with TRIS (hydroxymethyl) aminomethane (tris) buffer solution containing L-Aspartate, Malate Dehydrogenase (MDH) and Lactate Dehydrogenase (LDH). The reaction is initiated by the addition of NADH and measured kinetically. The rate of NADH consumption is directly proportional to the AST concentration in the patient sample.

INSTRUMENTS:

Instrument applications are available upon request.

COMPONENT COMPOSITION:

Component	Ingredients	Concentration in Tests
Reagent 1	tris Buffer pH 7.8	100 mmol/l
	L-Aspartate	300 mmol/l
	MDH	> 530 UI/l
	LDH	> 750 UI/l
	STABILISERS & PRESERVATIVES	—
Reagent 2	α -Ketoglutarate	75 mmol/l
	NADH	0.25 mmol/l
	PRESERVATIVES	—

REAGENT PREPARATION AND STABILITY:

Before use, mix reagent by gently inverting each bottle. If stored and handled properly, unopened components are stable until the expiry date stated on the label.

Monoreagent procedure: Add 1 volume of Reagent 2 to 5 Volumes of Reagent 1.

Working reagent is stable 5 weeks at 2-8°C.

Bireagent procedure: Liquid reagent 1 and 2 are ready for use.

Once open, components are stable until expiry date on label if store and handled properly.

TYPE OF SPECIMEN: 1

Use serum, free of haemolysis, heparin or EDTA plasma as specimen.

It is recommended to follow NCCLS procedures (or similar standardised conditions) regarding specimen handling. Specimen should be collected in an appropriate sample container, with proper specimen identification. Serum/plasma should be separated from cells within 48 hours after collection.

Stability: up to 5 days at 2-8°C³.

TEST PROCEDURE:

Materials required but not supplied:

Description	Catalog No.	Description	Catalog No.
General Chemistry Calibrator	GL983	Photometer	NA
General Chemistry Control Level 1	GL922	General Laboratory Equipment	NA
General Chemistry Control Level 2	GL932		

Assay procedure:

Wavelength: λ : 340 nm (Hg 334-Hg 365)

Temperature: 30°C or 37°C

Optical path: 1 cm light path.

CALCULATION:

340 nm: UI = AOD/min x 1746 334 nm: UI = AOD/min x 1780 365 nm: UI = AOD/min x 3235
(Conversion factor: Qty in μ Kat/l = Qty in UI x 0.0167).

EXPECTED VALUES:

IFCC reference method	30°C* UI	30°C* μ Kat/l	37°C* UI	37°C* μ Kat/l
Men	Up to 25	Up to 0.42	Up to 37	Up to 0.63
Women	Up to 21	Up to 0.36	Up to 31	Up to 0.53

*Calculated values *Temperature conversion factor from 37°C to 30°C: 0.67⁴
Each laboratory should establish its own reference range. AST results should always be reviewed with the patient's medical examination and history.

PERFORMANCE CHARACTERISTICS:

Performance results can vary with the instrument used. Data obtained in each individual laboratory may differ from these values.

Linearity:

Linear up to 412 UI (6.8 μ Kat/l).

For samples with a higher concentration, dilute 1:1 with 0.9% NaCl (9g/l) and re-assay. Multiply result by 2.

Interfering substances:

Bilirubin (mixed isomers): 12% interference up to 200 μ mol/l Bilirubin.
Haemolysis: 12% interference up to 1.25 g/l Haemoglobin.
Lipemia: Less than 10% interference up to 5 g/l Intralipid.

Sensitivity:

The Lowest Detectable Level was estimated at 3 UI (0.05 μ Kat/l).

Precision:

Within Run	Mean (UI)	SD	% CV	Between Run	Mean (UI)	SD	% CV
N = 20				N = 20			
Level 1	35.6	0.624	1.74	Level 1	40	1.09	2.74
Level 2	122.6	0.88	0.72	Level 2	133	2.94	2.22

Method Comparison:

Using 50 samples, a comparison, between this AST test (y) and another commercially available test (x), gave the following results:

$y = 1.014x - 0.219$	$r = 0.999$	Sample range: 9 to 355 UI.
----------------------	-------------	----------------------------

BIBLIOGRAPHY:

1. Bergmeyer HU et al. Clin Chem. 1971; 24:88 & Burles CA, Ashwood ER. Tetz Fund. Df Clin. Chem. 9th ed: 30-54 and 352-393.
2. Tietz NW, ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, Pa: WB Saunders, 1996: 76-77.
3. Fischbach F, Zawta B. Age-dependent Ref Limits of Several Etc in Plasma at Diff Measuring Temp. Klin Lab 1992; 38:585-591.
4. Zawta B, Klein G, Bialon W. Temperature Conversion in Clinical Enzymology Klin Lab 1994; 40:33-42.
5. IFCC Scientific Committee Clin. Chem. Biochem. 1990; 18:521-534.

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

Analytic reagents

Glucose oxidase	Glucose oxidation	<i>Aspergillus niger</i>	Detection of glucose in blood
Luciferase	Bioluminescence	Marine bacteria or firefly	Bioluminescent assays involving ATP
Peroxidase	Dye oxidation using H ₂ O ₂	Horseradish	Quantification of hormones and antibodies
Urease	Hydrolysis of urea to CO ₂ and NH ₃	Jack bean	Urea quantification in body fluids

Manipulative tools

Lysozyme	Hydrolysis of 1–4 glycosidic bonds	Hen egg white	Disruption of mucopeptide in bacterial cell walls
Nucleases	Hydrolysis of phosphodiester bonds	Various bacteria	Restriction enzymes used in genetic manipulation to cut DNA
DNA polymerases	DNA synthesis	<i>Thermus aquaticus</i>	DNA amplification used in the polymerase chain reaction

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
Composition of enzyme membranes available for analysers with a peroxide-sensitive electrode as the transducer.

Analyte	Enzyme	Reaction
Glucose	Glucose oxidase	$\beta\text{-D-glucose} + \text{O}_2 \rightarrow \text{gluconic acid} + \text{H}_2\text{O}_2$
Alcohol	Alcohol oxidase	$\text{Ethanol} + \text{O}_2 \rightarrow \text{acetaldehyde} + \text{H}_2\text{O}_2$
Lactic acid	Lactate oxidase	$\text{L-lactate} + \text{O}_2 \rightarrow \text{pyruvate} + \text{H}_2\text{O}_2$
Lactose	Galactose oxidase	$\text{Lactose} + \text{O}_2 \rightarrow \text{galactose dialdehyde derivative} + \text{H}_2\text{O}_2$

Biosensors = self-contained integrated devices which incorporate a biological recognition component (usually an immobilized enzyme) and an electrochemical detector (known as a transducer).

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



Horseradish Peroxidase Market: Segmentation

Grade


- *Grade I*
- *Grade II*
- *Grade III*

Application

- *Immunoassays*
 - *ELISA*
 - *Blot Test*
 - *Others*
- *Biosensors*

End User

- *Biopharmaceutical Companies*
- *Academic and Research Centers*
- *Others*


Substrate

- *ABTS*
(2,2'-Azinabis [3-ethylbenzothiazoline-6-sulfonic acid] -diammonium salt)
- *OPD*
(ophenylenediamine dihydrochloride)
- *TMB*
(3,3',5,5'-tetramethylbenzidine)
- *Others*



Region

- *North America*
- *Europe*
- *Asia Pacific*
- *Latin America*
- *Middle East & Africa*

Systemy s alkalickou fosfatázou

Plasmonic ELISA based on the controlled growth of silver nanoparticles


Alkaline phosphatase labeled antibody-based electrochemical biosensor for sensitive HE4 tumor marker detection

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
Assay components

- The components of **MEIA** include the following, all suspended in a specific buffer optimized for the assay:
 - Microparticle-Antibody Solid Phase:** Latex microparticles that are coated with antibody to bind the specific analyte being measured;
 - Antibody-Enzyme Conjugate:** Alkaline Phosphatase enzyme conjugated to antibody;
 - Enzyme Substrate:** Fluorescent 4-Methyl Umbelliferyl Phosphate (MUP) in solution that is available for a reaction with the enzyme on the antibody.

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



Image 11 Bayer System 240 bDNA (branched chain DNA) Analyzer.

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