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Oligo Calc: Oligonucleotide Properties Calculator

Enter Oligonucleotide Sequence Below
OD calculations are for single-stranded DNA or RNA

Nucleotide base codes

CAG TAT CCT GCC GAG TCT GT

Reverse Complement Strand(5' to 3') is:
ACA GAC TCG GCA GGA TAC TG

5' modification (if any) 3' modification (if any) Select molecule
50 nM Primer 1 Measured Absorbance at 260 nanometers ssDNA
50 mM Salt (Na⁺)

Calculate **Swap Strands** **BLAST** **mfold**

Physical Constants **Melting Temperature (T_M) Calculations**

Length: 20 Molecular Weight: 6084⁴ GC content: 55%
1 ml of a sol'n with an Absorbance of 1 at 260 nm
is 5.042 microMolar⁵ and contains 30.7 micrograms.


Thermodynamic Constants Conditions: 1 M NaCl at 25°C at pH 7.

RlnK 33.404 cal/(°K*mol) deltaH 161.9 Kcal/mol
deltaG 26.2 Kcal/mol deltaS 421.4 cal/(°K*mol)

Deprecated Hairpin/self dimerization calculations

5 (Minimum base pairs required for single primer self-dimerization)
4 (Minimum base pairs required for a hairpin)

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Format Conversion

- Combine FASTA
- EMBL to FASTA
- EMBL Feature Extractor
- EMBL Trans Extractor
- Filter DNA
- Filter Protein
- GenBank to FASTA
- GenBank Feature Extractor
- GenBank Trans Extractor
- One to Three
- Range Extractor DNA
- Range Extractor Protein
- Reverse Complement
- Split Codons
- Split FASTA
- Three to One
- Window Extractor DNA
- Window Extractor Protein

Sequence Analysis

Sequence Manipulation Suite:

Reverse Complement

Reverse Complement converts a DNA sequence into its reverse, complement, or reverse-complement counterpart. Reverse Complement converts a DNA sequence into its reverse, complement, or reverse-complement counterpart. Reverse Complement converts a DNA sequence into its reverse, complement, or reverse-complement counterpart. Reverse Complement converts a DNA sequence into its reverse, complement, or reverse-complement counterpart.

Paste the raw sequence or one or more FASTA sequences into the text area below. Input limit is 100,000,000 characters.

```
CGCCTGGAGAATATTTGGGATGAGACACCAGTGTATTTGCTCCAAGCAGCCTCTTTGACC
TAAACTTCC
AGGCAGGATCTTAAATGAAAAAGAGGTACAGGATGAGGAGAAAAACAAGAAATTTGGCCT
TTCTGTGGG
CCATCACTTGGGCAAGTCCAT
```

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Oligo Calc: Oligonucleotide Properties Calculator

Enter Oligonucleotide Sequence Below
OD calculations are for single-stranded DNA or RNA

Nucleotide base codes

CAG TAT CCT GCC GAG TCT GT

Reverse Complement Strand(5' to 3') is:
ACA GAC TCG GCA GGA TAC TG

5' modification (if any) 3' modification (if any) Select molecule
[] [] ssDNA

50 nM Primer
50 mM Salt (Na⁺) [1] Measured Absorbance at 260 nanometers

Calculate **Swap Strands** **BLAST** **mfold**

Physical Constants **Melting Temperature (T_M) Calculations**

Length: Molecular Weight: GC content: %
1 ml of a sol'n with an Absorbance of at 260 nm
is microMolar and contains micrograms.

Thermodynamic Constants Conditions: 1 M NaCl at 25°C at pH 7.

RlnK cal/(°K*⁴mol) deltaH Kcal/mol
deltaG Kcal/mol deltaS cal/(°K*³mol)

Deprecated Hairpin/self dimerization calculations

(Minimum base pairs required for single primer self-dimerization)
 (Minimum base pairs required for a hairpin)

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NCBI Resources How To Sign in to NCBI

Nucleotide Nucleotide Search Help

GenBank Send to: Change region shown Customize view

Homo sapiens NAD(P)H quinone dehydrogenase 1 (NQO1), transcript variant 1, mRNA

NCBI Reference Sequence: NM_000903.2
[FASTA](#) [Graphics](#)

Go to:

LOCUS NM_000903 2601 bp mRNA linear PRI 29-MAR-2018

DEFINITION Homo sapiens NAD(P)H quinone dehydrogenase 1 (NQO1), transcript variant 1, mRNA.

ACCESSION NM_000903

VERSION NM_000903.2

KEYWORDS RefSeq.

SOURCE Homo sapiens (human)

ORGANISM [Homo sapiens](#)
Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Primates; Haplorrhini; Catarrhini; Hominidae; Homo.

REFERENCE 1 (bases 1 to 2601)

AUTHORS Cheng X, Liu F, Liu H, Wang G and Hao H.

TITLE Enhanced glycometabolism as a mechanism of NQO1 potentiated growth of NSCLC revealed by metabolomic profiling

JOURNAL Biochem. Biophys. Res. Commun. 496 (1), 31-36 (2018)

PMID 29291405

Analyze this sequence Run BLAST **Pick Primers** Highlight Sequence Features Find in this Sequence Show in Genome Data Viewer

Articles about the NQO1 gene Redox modulation of NQO1. [PLoS One. 2018]

Enhanced glycometabolism as a mechanism of NQO1 po [Biochem Biophys Res Commun. 2018]

RNA-binding activity of TRIM25 is mediated by its PRY/SPRY domain and is requi [BMC Biol. 2017]

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Primer-BLAST *A tool for finding specific primers*

NCBI/ Primer-BLAST: Finding primers specific to your PCR template (using Primer3 and BLAST).

[Reset page](#) [Save search parameters](#) [Retrieve recent results](#) [Publication](#) [Tips for finding specific primers](#)

PCR Template

Enter accession, gi, or FASTA sequence (A refseq record is preferred) [Clear](#)

NM_000903.2

Or, upload FASTA file

Range

Forward primer [Clear](#)

Reverse primer

Primer Parameters

Use my own forward primer (5'->3' on plus strand) [Clear](#)

Use my own reverse primer (5'->3' on minus strand) [Clear](#)

PCR product size

Min	Max
<input type="text" value="70"/>	<input type="text" value="1000"/>

of primers to return

Primer melting temperatures (T_m)

Min	Opt	Max	Max T_m difference
<input type="text" value="57.0"/>	<input type="text" value="60.0"/>	<input type="text" value="63.0"/>	<input type="text" value="3"/>

Exon/intron selection

A refseq mRNA sequence as PCR template input is required for options in the section [Clear](#)

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Primer-BLAST *A tool for finding specific primers*

► NCBI/ Primer-BLAST: Finding primers specific to your PCR template (using Primer3 and BLAST).

[Reset page](#) [Save search parameters](#) [Retrieve recent results](#) [Publication](#) [Tips for finding specific primers](#)

PCR Template

Enter accession, gi, or FASTA sequence (A refseq record is preferred) [Clear](#)

Range

Forward primer	From	To	Clear
	196		
Reverse primer		1016	

Or, upload FASTA file

Primer Parameters

Use my own forward primer (5'→3' on plus strand) [Clear](#)

Use my own reverse primer (5'→3' on minus strand) [Clear](#)

PCR product size

Min	Max
200	500

of primers to return

Primer melting temperatures (T_m)

Min	Opt	Max	Max T _m difference
57.0	60.0	63.0	3

Exon/intron selection

A refseq mRNA sequence as PCR template input is required for options in the section

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Search mode Automatic

Database Refseq mRNA

Exclusion Exclude predicted Refseq transcripts (accession with XM, XR prefix) Exclude uncultured/environmental sample sequences

Organism 9606
Enter an organism name (or organism group name such as enterobacteriaceae, rodents), taxonomy id or select from the suggestion list as you type.
[Add more organisms](#)

Entrez query (optional)

Primer specificity stringency
Primer must have at least 2 total mismatches to unintended targets, including
at least 2 mismatches within the last 5 bps at the 3' end.
Ignore targets that have 6 or more mismatches to the primer.

Max target size 4000

Allow splice variants Allow primer to amplify mRNA splice variants (requires refseq mRNA sequence as PCR template input)

[Get Primers](#) Show results in a new window Use new graphic view

[Advanced parameters](#) Note: Parameter values that differ from the default are highlighted in yellow

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Search mode Automatic

Database Refseq mRNA

Exclusion Exclude predicted Refseq transcripts (accession with XM, XR prefix) Exclude uncultured/environmental sample sequences

Organism 9606
Enter an organism name (or organism group name such as enterobacteriaceae, rodents), taxonomy id or select from the suggestion list as you type.

Entrez query (optional)

Primer specificity stringency
Primer must have at least 2 total mismatches to unintended targets, including
at least 2 mismatches within the last 5 bps at the 3' end.
Ignore targets that have 6 or more mismatches to the primer.

Max target size 4000

Allow splice variants Allow primer to amplify mRNA splice variants (requires refseq mRNA sequence as PCR template input)


Get Primers Show results in a new window Use new graphic view

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Input PCR template [NM_000903.2](#) Homo sapiens NAD(P)H quinone dehydrogenase 1 (NQO1), transcript variant 1, mRNA
Range 196 - 1016
Specificity of primers Primer pairs are specific to input template as no other targets were found in selected database: Refseq mRNA (Organism limited to Homo sapiens)
Other reports [Search Summary](#)

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Products on allowed transcript variants

>[NM_001025434.1](#) Homo sapiens NAD(P)H quinone dehydrogenase 1 (NQO1), transcript variant 3, mRNA

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