

STAPHYLOCOCCI,
STREPTOCOCCI

PRACTICAL PART

1. IDENTIFICATION OF STAPHYLOCOCCI


Free plasmacroagulase

- **Principle:** plasmacroagulase enzyme converts plasma soluble fibrinogen into insoluble fibrin
- **Tools:** test-tube with 0.5 mL of rabbit plasma, fresh staphylococcal culture, loop
- **Procedure:** Inoculate several staphylococcal colonies into the test tube with plasma, incubate 24 hours at 37°C
- **Evaluation:** Observe clot (coagulation) formation after 1 and 24 hours




Bound coagulase

- **Principle:** bound coagulase (Clumping factor) is part of the surface of coagulase-positive staphylococci – It leads to clumping of staphylococci and their protection against phagocytosis
- **Tools:** staphylococcal culture, plasma, microscopic slide, loops
- **Procedure:** drop rabbit plasma on the slide, add staphylococcal culture. Mix it properly.
- **Evaluation:** observe agglutination in case of *S. aureus*



Catalase test




- **Principle:** using for differentiation between staphylococci from streptococci
- Staphylococci produce catalase enzyme (disintegration of hydrogen peroxide – bubbles formation). Typical feature of aerobic or facultative anaerobic bacteria.
- **Tools:** culture of staphylococci and streptococci, hydrogen peroxide, microscopic slide, loops, disinfection
- **Procedure:** Add few colonies of bacterial culture on the slide than drop hydrogen peroxide onto culture. Observe bubbles formation.
- **Evaluation:** observe positivity (bubbles formation) in staphylococci, negativity in streptococci

PRACTICAL PART

2. IDENTIFICATION OF STREPTOCOCCI

Optochin test

- **Principle:** *Streptococcus pneumoniae* is differentiated from other viridans streptococci on the basis of its susceptibility to antibiotic - optochin (Ethyl Hydrocuprein)



Optochin test – practice


- **Tools:** *Streptococcus pneumoniae* and viridans streptococci cultures, blood agar, optochin discs, loops, markers, disinfection
- **Procedure:** inoculate several colonies of *Str. pneumoniae* on the first half of blood agar, viridans streptococci on the second half. Cultures of both strains must not overlapped!!!
- Put optochin disc on each strains and incubate 18-24 hours at 37°C, in 5% CO₂
- **Evaluation:** measure the zone of inhibition, in case of more than 14 mm report as susceptible

PRACTICAL PART

3. IDENTIFICATION OF BETA-HEMOLYTIC STREPTOCOCCI – CAMP TEST

CAMP test

- **Principle:** group B streptococci (*Streptococcus agalactiae*) produce extracellular protein (CAMP factor), which is able potentiate beta-hemolysin of *Staphylococcus aureus*



CAMP test – practice

- **Tools:** bacterial culture of tested strain, blood agar, *St. aureus* strain, marker, loops, disinfection
- **Procedure:** inoculate the line of tested streptococcal strain on blood agar and perpendicularly to this line inoculate the line of *St. aureus*
- Incubation 24 hours at 37°C, aerobically

PRACTICAL PART

4. DIFFERENTIATION OF ENTEROCOCCI FROM STREPTOCOCCI

Bile - esculin test

Principle:

- Differentiation of enterococci from other streptococci using cultivation in highly selective-diagnostic bile – esculin medium
- Bile inhibits the growth of streptococci
- Enterococci hydrolyzes esculin and medium color turns black

Bile - esculin test

Tools: bacterial culture of tested strain, bile-esculin medium, marker, loops


Procedure:

- Pick tested bacterial culture using the loop and inoculate via performing the puncture through the medium to the bottom of tube.
- Incubation 24 hours at 37°C, aerobically

Evaluation:

Positive reaction – **blackening of medium**

Negative reaction – color of medium stays unch




PRACTICAL PART

5. IDENTIFICATION OF STREPTOCOCCI – EVALUATION

CAMP test

- **Evaluation:** positive test is shown as enhancement of beta-hemolysis of tested strain (butterfly wings like) in case of group B streptococci



Optochin test

- **Evaluation:** Measure the zones of inhibition around antibiotics after the incubation. Notice the diameter and interpret result of identification according to optochin susceptibility



Inhibition zone \geq 14 mm

Report as susceptible

Streptococcus pneumoniae

Inhibition zone $<$ 14 mm

Report as resistant

viridans streptococci

PRACTICAL PART

6. IDENTIFICATION OF VIRIDANS

STREPTOCOCCI – BILE SOLUBILITY TEST

Confirmation of identification of *Streptococcus pneumoniae* – bile solubility test

Principle:

- Deoxycholate (bile) activates autolysin in cell wall of *Str. pneumoniae*
- *Str. pneumoniae* colonies are dissolved in presence of deoxycholate
- Difference – viridans streptococci colonies show no change in presence of deoxycholate (absence of autolysin in cell wall)

Bile solubility test - practice

- **Tools:** *Str. pneumoniae* and viridans streptococci colonies, 10% solution of sodium deoxycholate
- **Procedure:** Drop deoxycholate on tested colony (colonies). After 15 min observe and report if colony was dissolved or not.
- **Evaluation:**

Dissolved colonies – report as *Str. pneumoniae*

Colonies are intact – viridans streptococci

PRACTICAL PART

7. LATEX AGGLUTINATION OF BETA- HEMOLYTIC STREPTOCOCCI

Detection of group antigen of beta-hemolytic streptococci

Principle:

Typing of beta-hemolytic streptococci according to specific type of polysaccharide C in streptococcal cell wall (Lancefield antigen). Commercial sets – known types of antibodies (against different types of Lancefield antigens) bound on latex particles.

Antigen presence will appear as agglutination with specific antibody

Tools: Group specific antisera (antibodies on latex particles), cell wall extract of streptococcal strain, plates, droppers, sticks, disinfection


Latex agglutination of beta-hemolytic streptococci

- **Procedure:** transfer one drop of extract (prepared from 5 colonies of tested strain in advance) to each of the five black circles labeled A, B, C, F, G.
- Shake the bottles with latex particles anti-A, B, C, F and G and add one drop close to the drop of extract according to the label. Mix the contents in each circle in turn with a stick and spread to cover area of the circle
- Rock the card gently for a maximum of 1 minute and read the reaction.

Latex agglutination of streptococci

Evaluation:

- Observe the formation of agglutination during the rotation



- <http://www.pro-lab.com/products-strep.php?country=DK>
- <https://www.youtube.com/watch?v=xicBT3FaSM0>

PRACTICAL PART

8. EVALUATION OF BILE-ESCULIN TEST

Bile - esculin test

After 24 hours incubation at 37°C, aerobically evaluate the colour change of medium

Evaluation:

Positive reaction – **blackening of medium**

Negative reaction – color of medium stays unchanged

